
SMART 2.0 Operation Manual

Data : 2023.05.24

Directory

| | |
|---------------------------------|---|
| 1. Software installation | 3 |
| 1.1 Operating Environment | 3 |
| 1.2 Hardware Requirements | 3 |
| 1.3 Installation steps | 3 |
| 2. Data input format | 3 |
| 2.1 Data Processing | 3 |
| 2.2 Data Types | 4 |
| 3. Steps for using SMART | 5 |
| 3.1 Create a new project. | 5 |
| 3.2 Input raw data | 6 |
| 3.3 Profile processing | 6 |
| 3.4 Resolve Profile | 7 |
| 3.5 Calculate LR | 8 |
| 4. View results | 8 |
| 4.1 View results in smart | 8 |
| 4.2 View results in pdf | 9 |

1. Software installation

1.1 Operating Environment

Operating System: Windows 7, Windows 10, Windows 11

System digits: 64-bit

1.2 Hardware Requirements

CPU: X86-based

Memory: 16GB or more

Disk space: 20G or more

1.3 Installation steps

- a) Download smart_2.0.zip.
- b) Copy smart_2.0.zip to the installation path, such as: C:\, and unzip it.
- c) Open C:\smart_2.0\GenerateCmd\
- d) Double-click the file 'get-C2V.exe' to complete the installation.
- e) Open C:\smart_2.0, double-click the Smart icon to launch it.

2. Data input format

2.1 Data Processing

SMART can directly process the offline data from genetic analyzers, or process the electrophoresis graphs in fluorescence mapping software. During the process, OL peaks, penetration peaks, sharp peaks, and artificially identified insertion peaks need to be removed. Note that SMART has already modeled the shadow peaks at the minus one position (backward) and plus one position (forward) of the allele, so there is no need to delete the forward and backward shadow peaks. However, other types of shadow peaks, such as those at the minus two position or plus two position, need to be deleted. For some specific loci (e.g., SE33 locus), shadow peaks that are two base positions shorter may also need to be deleted.

After the graph processing is completed, the graph data is exported as a text file (TXT format) in the fluorescence mapping software (taking GenemapperID-X

software as an example).

2.2 Data Types

The input data for SMART includes two types:

(1) Data files of mixed DNA samples, consisting of five main parts: sample name, Marker name, allele name, size of the allele, and corresponding peak height of the allele.

| Sample Name | Marker | Allele 1 | Allele 2 | Allele 3 | Allele 4 | Allele 5 | Allele 6 | Allele 7 | Allele 8 | Allele 9 | Size 1 | Size 2 | Size 3 | Size 4 | Size 5 | Size 6 | Size 7 | Size 8 | Size 9 | Height 1 | Height 2 | Height 3 | Height 4 | Height 5 | Height 6 | Height 7 | Height 8 | Height 9 | | | |
|-------------|---------|----------|----------|----------|----------|----------|----------|----------|----------|----------|--------|--------|--------|--------|--------|--------|--------|--------|--------|----------|----------|----------|----------|----------|----------|----------|----------|----------|--|--|--|
| test | D8S1179 | 9 | 10 | 12 | 13 | 14 | 15 | 16 | | | 125.64 | 129.66 | 137.82 | 141.99 | 146.18 | 150.3 | 154.44 | | | 171 | 3665 | 119 | 1953 | 114 | 306 | 4086 | | | | | |
| test | D21S11 | 29 | 30 | 31 | 32 | | | | | | 203.01 | 207.01 | 210.93 | 215.01 | | | | | | 448 | 639 | 92 | 782 | | | | | | | | |
| test | D7S820 | 7 | 8 | 10 | 11 | 12 | | | | | 258.32 | 262.29 | 270.35 | 274.29 | 278.37 | | | | 69 | 1500 | 551 | 283 | 1918 | | | | | | | | |
| test | CSF1PO | 9 | 11 | 12 | | | | | | | 315.02 | 323.43 | 327.46 | | | | | | 253 | 169 | 1567 | | | | | | | | | | |
| test | D3S1358 | 14 | 15 | 16 | 17 | 18 | | | | | 118.35 | 122.27 | 126.24 | 130.26 | 134.43 | | | | 432 | 5393 | 516 | 3926 | 572 | | | | | | | | |
| test | TH01 | 6 | 7 | 8 | 9 | 9.3 | | | | | 168.83 | 172.9 | 176.98 | 181.1 | 184.17 | | | | 103 | 6507 | 100 | 2657 | 6223 | | | | | | | | |
| test | D13S317 | 7 | 8 | 9 | 11 | 12 | 13 | | | | 211.15 | 215.13 | 219.2 | 227.26 | 231.31 | 235.36 | | | 54 | 2074 | 269 | 51 | 629 | 1236 | | | | | | | |
| test | D16S539 | 9 | 10 | 11 | | | | | | | 267.46 | 271.51 | 275.45 | | | | | | 730 | 2829 | 3416 | | | | | | | | | | |
| test | D2S1338 | 17 | 18 | 19 | 20 | 23 | 24 | | | | 313.16 | 317.44 | 321.6 | 325.63 | 337.6 | 341.61 | | | 85 | 1329 | 1444 | 113 | 74 | 173 | 2232 | | | | | | |
| test | D19S433 | 12 | 12.2 | 13 | 13.2 | 14 | 14.2 | 15 | | | 112.51 | 114.58 | 116.52 | 118.55 | 120.49 | 122.56 | 124.45 | | 461 | 186 | 7057 | 2399 | 6183 | 213 | 2232 | | | | | | |
| test | vWA | 14 | 15 | 16 | 17 | | | | | | 164.35 | 168.19 | 172.25 | 176.34 | | | | | 1394 | 104 | 1783 | 296 | | | | | | | | | |
| test | TPOX | 7 | 8 | 9 | 11 | | | | | | 225.02 | 229.06 | 233.11 | 241.17 | | | | | 53 | 2270 | 1986 | 338 | | | | | | | | | |
| test | D18S51 | 12 | 13 | 14 | 15 | 16 | 17 | | | | 281.28 | 285.34 | 289.29 | 293.36 | 297.32 | 301.25 | | | 72 | 1476 | 3696 | 138 | 1672 | 3754 | | | | | | | |
| test | AMEL | X | Y | | | | | | | | 104.92 | 110.92 | | | | | | | 1962 | 2424 | | | | | | | | | | | |
| test | D5S818 | 11 | 12 | 13 | | | | | | | 148.92 | 153.06 | 157.21 | | | | | | 139 | 1150 | 1273 | | | | | | | | | | |
| test | FGA | 20 | 21 | 22 | 23 | 24 | | | | | 224.91 | 228.95 | 232.99 | 237.06 | 241.17 | | | | 672 | 163 | 2999 | 1468 | 3740 | | | | | | | | |

Figure 2.1 Data files of mixed DNA samples

(2) Data files of single-person DNA samples, which are STR typing files for a single individual. These are generally STR typing files for suspects or typing files for known contributors of mixed samples. The file format is the same as the mixed DNA file format.

| Marker | Allele 1 | Allele 2 | Allele 3 | Allele 4 | Allele 5 | Allele 6 | Allele 7 | Allele 8 | Allele 9 | Size 1 | Size 2 | Size 3 | Size 4 | Size 5 | Size 6 | Size 7 | Size 8 | Size 9 | Height 1 | Height 2 | Height 3 | Height 4 | Height 5 | Height 6 | Height 7 | Height 8 | Height 9 |
|---------|----------|----------|----------|----------|----------|----------|----------|----------|----------|--------|--------|--------|--------|--------|--------|--------|--------|--------|----------|----------|----------|----------|----------|----------|----------|----------|----------|
| D8S1179 | 9 | 10 | | | | | | | | | 125.64 | 129.66 | | | | | | | 171 | 3665 | | | | | | | |
| D21S11 | 29 | | | | | | | | | | 203.01 | | | | | | | | 448 | | | | | | | | |
| D7S820 | 7 | 8 | | | | | | | | | 258.32 | 262.29 | | | | | | | 69 | 1500 | | | | | | | |
| CSF1PO | 9 | | | | | | | | | | 315.02 | | | | | | | | 253 | | | | | | | | |
| D3S1358 | 14 | 15 | | | | | | | | | 118.35 | 122.27 | | | | | | | 432 | 5393 | | | | | | | |
| TH01 | 6 | 7 | | | | | | | | | 168.83 | 172.9 | | | | | | | 103 | 6507 | | | | | | | |
| D13S317 | 7 | 8 | | | | | | | | | 211.15 | 215.13 | | | | | | | 54 | 2074 | | | | | | | |
| D16S539 | 9 | | | | | | | | | | 267.46 | | | | | | | | 730 | | | | | | | | |
| D2S1338 | 17 | | | | | | | | | | 313.16 | 317.44 | | | | | | | 85 | 1329 | | | | | | | |
| D19S433 | 12 | 18 | | | | | | | | | 112.51 | | | | | | | | 461 | | | | | | | | |
| vWA | 14 | 15 | | | | | | | | | 164.35 | 168.19 | | | | | | | 1394 | 104 | | | | | | | |
| TPOX | 7 | 8 | | | | | | | | | 225.02 | 229.06 | | | | | | | 53 | 2270 | | | | | | | |
| D18S51 | 12 | | | | | | | | | | 281.28 | | | | | | | | 72 | | | | | | | | |
| AMEL | X | Y | | | | | | | | | 104.92 | 110.92 | | | | | | | 1962 | 2424 | | | | | | | |
| D5S818 | 11 | | | | | | | | | | 148.92 | | | | | | | | 139 | | | | | | | | |
| FGA | 20 | 21 | | | | | | | | | 224.91 | 228.95 | | | | | | | 672 | 163 | | | | | | | |

Figure 2.2 Sample of single DNA sample data file

3. Steps for using SMART

3.1 Create a new project.




Figure 3.1 Main interface of SMART

Click on “Create a new project” to create a new task.





Figure 3.2 Create a new project


Enter Project ID and Sample ID. Then, click on “Create”.

3.2 Input raw data




Click on “Open hid/fsa folder” and choose a folder that contains the fluorescence data files to be analyzed (files with the ".hid" suffix). This folder should contain at least one file named "Ladder". Then click on “Analyze Profile”.

3.3 Profile processing



Process the profile, delete spurious peaks except allelic peaks and stutter peaks, and click save after processing.


3.4 Resolve Profile



1. Click on “Resolve Profile”
2. Click on “Open”
3. Select the mixed file that has just been processed. If there are parallel amplification files, multiple parallel amplification files can be selected and inputted simultaneously.
4. Click on “Confirm”

If there are known contributors in the mixed samples, they can be added by clicking on “Input profile of known contributor”. Then choose “Kit type” and input “num of contributors”. Finally click on “Run”.

3.5 Calculate LR




After running "Resolve Profile", in order to calculate LR, click "Calculate LR", add the POI's genotype, and then click "Run".

4. View results

4.1 View results in smart

Click on "Results".



On the left side of this interface, there are the results of the resolve, which allows

you to view the genotyping results of different contributors (Contributor 1~N) at different loci that exceed the specified threshold.

On the upper right side of this interface, there are the results of the Likelihood Ratio (LR) calculations. Each row represents the LR calculation result for a particular Person of Interest (POI), including information such as the LR value and the sample name.

4.2 View results in pdf

The detailed calculation results of resolution are located under
*/smart_results/ProjectID_SampleID_time/Reports/Resolve.pdf.

The detailed calculation results of LR are located under
*/smart_results/ProjectID_SampleID_time/Reports/ LR_POI1.pdf.